ARGON MEDICAL DEVICES, INC.
 1445 Flat Creek Road Athens, Texas 75751
 USA Tel: +1 (903) 675 9321
 Tel: +1 (800) 927 4669
 www.argonmedical.com

EMERGO EUROPE
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands
 +31 70 345 8570

Argon Medical Devices UK Ltd
 Eastgate Business Centre
 Eastern Avenue
 Burton-on-Trent
 Staffordshire
 DE13 0AT

IFU8000/0321A

<https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Portuguese	O glossário dos símbolos está localizado eletronicamente em www.argonmedical.com/symbols
Czech	Slovníček symbolů je elektronicky umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse electrónicamente en www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols .
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Croatian	Pojmovnik simbola nalazi se u elektroničkom obliku na www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Lithuanian	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Russian	Электронный глоссарий символов расположен по адресу www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols
Swedish	Symbolforklaringen hittas elektroniskt på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır
Arabic	www.argonmedical.com/symbols يتوفر مسرد الرموز إلكترونياً على

Intended Use/Purpose

The Polyethylene (PE) Drainage Catheter product family is intended for use in general drainage.

- TCD Drainage Catheters and Sets are for drainage of cysts, abscesses, hematomas, plural exudates, ascites, gall bladders, nephrostomies, pleuracentesis and paracentesis.
- TCD Drainage Catheters with Safety String Lock System are for drainage of cysts, gall bladders and nephrostomies.
- Choice Lock™ TCD Drainage Sets are for single step drainage of cysts, abscesses, hematomas, plueral exudates, ascites, gall bladders and nephrostomies.
- Choice Lock™ TCD Drainage Sets with Safety String Lock System are for single step drainage of cysts, gall bladders and hematomas. PE Biliary Catheters and Sets are for biliary drainage.
- PE Nephrostomy Catheters (pigtail) are for drainage cysts, gall bladders and nephrostomies.
- Nephrostomy Drainage Catheter Sets are for external urine drainage from the renal pelvis.

Device Description

- These devices are considered a long-term surgical invasive device permanently implanted for more than 30 days.
- These devices are not intended for direct contact with the heart or central circulatory system or central nervous system.
- TCD drainage catheters are devices made of radiopaque polyethylene material, designed to be used in a variety of drainage applications and come with or without a safety string lock system. The Single step version comes with the Choice Lock® needle which allows for simple singlehanded snap lock and release of the stylet.
- PE Biliary drainage catheters are introduced percutaneously for temporary or prolonged drainage of the biliary tract for splinting of the bile duct during healing or preventing stricture of the bile duct.
- All drainage catheters are sterilized with a Sterility Assurance Level (SAL) of 10⁻⁶ after exposure to the ethylene oxide gas sterilization process in conformance to consensus standard ISO 11135- 1:2014, Sterilization of health care products – Ethylene Oxide – Part 1: Requirements for development, validation and routine control of sterilization process for medical devices.

Indication for Use

The All-Purpose Drainage Catheters and Sets are single-use devices intended for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary and abscess), either by direct stick or Seldinger access technique. Surgical drains are used to decompress or drain either fluid or air from the area of surgery to prevent the accumulation of fluid or air.

Contra-indications for Use

There are no absolute contraindications with use of this product. However, the assessment of the patient's overall physical state by the physician is required prior to the procedure.

Complications/Side Effects

- Leakage
- Tract Infection
- Allergic reaction
- Temporary pain or discomfort
- Infection
- Bleeding
- Nontarget Puncture
- Hemorrhage
- Sepsis
- hepatocolic fistula
- Pseudoaneurysm
- Adjacent vessel laceration

***Drain selection and placement is important to ensure that possible complications are minimized. Drains may become occluded or clogged, resulting in retained fluid. Patency must be maintained.

Side Effects

There are no known side effects.

Multipurpose Drainage Catheters and Sets

Warnings

- These devices were designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Do not use if package appear open or if the expiry date has been exceeded.
- The product must only be used by qualified personnel who are familiar with the technique.
- During insertion avoid contact with bone, cartilage and scar tissue which can damage the catheter tip.

Preparation for Devices with Catheters or Pigtailed

- Advance the straightener carefully over the curvature of the catheter, while straightening the curvature with your fingers.
- Gently pull the thread to lock the pigtail in order to avoid unintentional looping of the thread.
- Introduce the chosen stiffener (if plastic: activate the hydrophilic coating with saline) and fasten the stiffener onto the catheter (Luer Lock).
- Remove the straightener
- Activate the coating by saline.

Procedure – Biliary Drainage Catheter

Placement

- Advance the catheter over the guidewire that has been introduced into the patient using standard technique.
- When the catheter is fully loaded onto the guidewire remove the straightener.
- Verify correct position of the catheter using fluoroscopy.
- Remove the guidewire.
- Connect the catheter to a drainage bag using a suitable connecting tube.
- We recommend bandaging the catheter in a straight line as it exits the skin. Any curvature required can be applied to the connecting tube.

Procedure – TCD Drainage Catheter

Placement Using Catheter

- Advance the catheter over the guidewire that has been introduced into the patient using standard technique.
- When the catheter is fully loaded onto the guidewire remove the straightener.
- Ensure that the complete length of the pigtail is in the cavity.
- Loosen the adapter and gently pull the thread turning the catheter clockwise to ensure correct position and to complete curvature.
- Fix the thread tightly by winding it around the slot in the adapter and finish up with pressing the clip onto the slot and twist off the handle from the clip.
- Verify correct position of the catheter using fluoroscopy.
- Remove the guidewire
- Connect the catheter to a drainage bag using a suitable connecting tube.
- It is recommended affixing the catheter in a straight line and that any curvature be applied to the connecting tube.
- It is recommended that the physician follow the hospital's standard of care procedures for drainage catheters.

Removal of Catheter

- Disconnect the drainage bag connector tube from catheter.
- Loosen the female/male adapter from the standard luer hub at the catheter to deactivate the valve.
- Remove the clip and unwind the suture. Check that both threads are loose and cut one thread in order to loosen the pigtail.
- Pull the catheter out gently. If access is to be maintained, a straight floppy tip guidewire passed through the catheter will facilitate removal while maintaining access.

Caution for Locking Catheters

- If the catheter is to be removed on another ward, we recommend that these guidance notes accompany the patient's case notes to ensure that the relevant personnel are aware of the presence of a locking catheter. It is also advisable to inform the patient.

Procedure – Nephrostomy Drainage Catheter

Placement Using Catheter

- Advance the catheter over the guidewire that has been introduced into the patient using standard technique.
- When the catheter is fully loaded onto the guidewire

remove the straightener.

- Advance the catheter to the drainage site
- Verify correct position of the catheter using fluoroscopy.
- Remove the guidewire.
- Connect the catheter to a drainage bag using a suitable connecting tube.
- It is recommended affixing the catheter in a straight line and that any curvature be applied to the connecting tube.
- It is recommended that the physician follow the hospital's standard of care procedures for drainage catheters.

Removal of Catheter

- Disconnect the drainage bag connector tube from catheter.
- Remove the clip and unwind the suture. Check that both threads are loose and cut one thread in order to loosen the pigtail.
- Pull the catheter out gently. If access is to be maintained, a straight floppy tip guidewire passed through the catheter will facilitate removal while maintaining access.

Duration

Up to 8 weeks.

Storage

Store at controlled room temperature. Avoid extended exposure to light.

Disposal

After use, this product may be a potential biohazard. Dispose in accordance with applicable laws and regulations.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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