ARGON

MEDICAL DEVICES

ARGON MEDICAL DEVICES, INC.
 1445 Flat Creek Road, Athens, Texas 75751 USA
 Tel: +1 800-927-4669;
 Tel: +1 903-675-9321
www.argonmedical.com

EMERGO EUROPE
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands
 +31 70 345 8570

Argon Medical Devices UK Ltd

Eastgate Business Centre
 Eastern Avenue
 Burton-on-Trent
 Staffordshire
 DE13 0AT

99-930158010/0221A

 <https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Bulgaria	Речникът на символите е наличен в електронен вид на адрес www.argonmedical.com/symbols
Czech	Slovníček symbolů je elektronicky umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse electrónicamente en www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols .
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Lithuanian	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Latvian	Simbolu glosārijs ir atrodams elektroniskā veidā vietnē www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Slovenian	Glosár simbolov je v elektronickej podobi k dispoziciji na stránkach www.argonmedical.com/symbols
Swedish	Symbolforklaringen hittas elektroniskt på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır

INTENDED USE/PURPOSE

This product is used in angiographic or radiologic procedures in controlling fluids/waste and the injection of contrast media into the patients.

DESCRIPTION

Patients Monitoring Products are stand-alone components, stock kits, or convenience kits built to user specifications. Products addressed by these instructions for use include: monitoring sets, pressure transducers, patient mount straps, continuous flush devices, stopcocks, manifolds, arterial extensions, IV accessories, needle-less sample site, hyperbaric pass through and infusion kits, KeepSafe™ Closed Sampling System options, contrast line and Argo Bag.

Contents of these sets have been sterilized with ethylene oxide gas. They are only for single patients, single-use applications.

INDICATIONS FOR USE: Patient Monitoring products are indicated for physiological pressure measurement, invasive pressure monitoring, catheterization procedures, infusion of fluids, and/or obtaining blood or other fluid samples. Contrast Line are indicated to conduit tubing to deliver fluids from an electronic or manual injector. Connector tubing to provide a sterile fluid pathway between two devices.

CONTRAINDICATIONS: There are no known contraindications for these products.

Possible complications:

There no known complications related to these devices when used as intended.

WARNINGS

- Carefully read these instructions before using these products. When using sets incorporating other manufacturer's components, also read their current instruction for use.
- This device is designed and intended for ONE-TIME USE ONLY. Do not resterilize and/or reuse. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. This device is extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.
- Do not use Patient Monitoring Products with parenteral nutritional solutions containing lipids or nitroglycerin. These substances are known to degrade some medical plastics resulting in cracks and/or leaks. Other drugs or solutions may have similar effects. Before using drugs or solutions with which you are unfamiliar, please consult the drug manufacturer's IFU for plastics compatibility information.
- Remove all air from the IV solution bag, all monitoring components and any reservoir syringes during setup. This reduces the possibility of air entering the patient's vascular system from the kit or IV bag if it runs out of fluid.
- Inspect for damage, proper assembly, and loose connections prior to use. Check for leaks before and during use. Leaks can result in the loss of sterility, fluid or blood loss, and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the leaking component. Do not over-tighten connections; all connections should be finger-tight. Over-tightened connections may result in cracks or leaks. Wet connections promote over-tightening by lubricating the ports.
- Do not use kits or components for intracranial pressure monitoring unless they are labeled as ICP product.
- For intraventricular, percutaneous, or intracranial pressure monitoring applications, follow the correct procedures for your institution for controlling drainage of intraventricular cerebrospinal fluid to prevent contamination of the system.
- Do not use a continuous flush device for intracranial pressure monitoring.
- Do not use a continuous flush device on neonatal or pediatric patients without using an IV flow controller.
- Only use the 30mL/hr continuous flush device (yellow band) with an infusion pump or other flow system between the continuous flush device and the fluid source. Do not use the 30mL/hr continuous flush device as a flow controller.
- When using a continuous flush device, always use a backup flow indicator such as a standard IV drip chamber, use this to check for a flush device that may be occluded or stuck in the fast-flush mode. Check this at set-up and periodically during use. Also, monitor the bag cuff pressure periodically to maintain a constant flow.
- Fast-flushing may not be appropriate for neonatal or pediatric monitoring or for other applications where patient volume control is necessary.
- Fluids with a viscosity greater than saline may reduce the continuous flush rate of the continuous flush device or cause it to occlude.
- Do not exceed the pressures below when using the following components:
 - Poly Ethylene Monitoring Lines and Rotating Adapters: 25,000 mmHg (500 psi)
 - PVC Monitoring Lines: 10,000 mmHg (200 psi)
 - Hyperbaric Pass Through Kits: 3ATA chamber pressure
 - Needle-less Sample Site: 1500 mmHg (29 psi)
 - Stopcocks: 775 mmHg (15 psi)
 - All other components: 500 mmHg (10 psi)

- For products containing the PHT symbol on the label: This product contains Di(2-ethylhexyl)phthalate (DEHP) which has been shown to cause reproductive harm in male neonates, pregnant women carrying male fetuses, and peripubertal males. The following procedures have been identified as posing the greatest risk for DEHP exposure: exchange transfusion in neonates, total parenteral nutrition (TPN) in neonates (with lipids in polyvinylchloride (PVC) bag), multiple procedures in sick neonates (high cumulative exposure), heart transplantation or coronary artery bypass graft surgery (aggregate dose), and massive infusion of blood into trauma patient. It is recommended that

DEHP-free medical products be considered when these procedures are to be performed on male neonates, pregnant women who are carrying male fetuses, and peripubertal males.

CAUTIONS

1. When using the needle-less sample site NOT necessarily all ports are closed when handle is turned 45° from ports.
2. When using the needle-less sample site malfunctioning of the valve might occur if non-standard luer is used.
3. Do NOT use needles to access the valve port of the needle-less sample site.
4. When using the needle-less sample site replace the site every 24 hours or per CDC guidelines and/or hospital protocol.

PRECAUTIONS

1. Do not use damaged products or products from packages that have tears, punctures, or broken seals.
2. Use proper aseptic techniques while handling these products.
3. Eliminate all air bubbles and align all stopcock handles exactly on the center of the port for the best possible response. The presence of air in the system will result in poor transmission of the patient's pressure.
4. These products are only for use by trained personnel only.
5. Maintain the drip chamber at least half full to minimize the risk of air entering the system while fastflushing.
6. Use the correct accessory cable to connect the transducer to the monitor.
7. Be sure that the transducer and cable connectors are dry, aligned, tight, and locked with the locking ring. Wet or improper connections may cause inaccurate readings.
8. Avoid kinking or cutting the electrical cables. Do not roll carts over cables.

TRANSDUCER SPECIFICATIONS

Defibrillation: Withstands 5 discharges of 400 joules in 5 minutes across 50 ohm load.

Specifications	CDX III & CDXpress	ArgoTrans
Pressure range	-50 to +300mmHg	-30 to +300mmHg
Overpressure	-500 to 5000 mmHg	-500 to 5000 mmHg
Operating temperature range	+15°C to +40°C	+10°C to +40°C
Storage temperature range	-30°C to +60°C	-25°C to +70°C
Eight hour drift	1mmHg after 10 min. warm-up	2mmHg after 10 min. warm-up

PROCEDURES FOR USE

Monitoring Sets

When using a monitoring set observe the following:

1. Connect the transducer to the accessory cable for the monitor being used. Align and firmly join the connectors together. Finger-tighten the locking ring; then connect the cable to the monitor. For greater accuracy, allow a minimum of five minutes warm-up time after connecting the transducer before zeroing.
2. To help maintain catheter patency, the IV solution may be heparinized per hospital protocol. For kits with an IV set, close the roller clamp and insert the IV spike into the solution bag. Without pressurizing the solution bag, fill drip chamber half full by gently squeezing the drip chamber. For kits without an IV set, attach IV solution source to the Luer lock connector at the end of the flush device. Elevate the bag to prime the system.
3. Mount the transducer according to hospital protocol.
4. The continuous flush device is designed to continuously infuse pressurized solutions. The 3 mL/hr device is for adult continuous infusion with an inlet pressure of 300 mmHg. The 30 mL/hr neonatal/pediatric models (yellow band) are only for use with an infusion pump or other flow controller. It delivers a nominal flow of 30 mL/hr (yellow band) with an inlet pressure of 300 mmHg.

When using the Continuous Flush Device, observe the following:

- Connect the infusion fluid source to the Luer lock connector at the end of the continuous flush device.
 - To prime the system or supply a fast flow of infusion fluid to the patient, either pull the flush tab or squeeze the wings on the continuous flush device. Release the pull tab or wings to return to the continuous flush mode.
5. After completely priming and removing all air from the system, place the solution bag in a pressure cuff. Position the transducer at the patient's normal zero reference level and pressurize the infusion cuff to 300 mmHg. If using an infusion pump or other method of controlling infusion, set up the infusion control system per manufacturer IFU.
 6. Connect the monitoring set to the patient's catheter without allowing air into the system. Flush residual blood and remaining air from the catheter. Turn on the in-line stopcock(s) to the patient.
 7. To zero balance the system (make sure the transducer is lower than the zero port cap), turn the stopcock in the patient line OFF to the patient. Open the zero port cap about ¼ turn to atmospheric pressure. Calibrate the monitor to 0 mmHg on the monitor display. Close the zero port cap and turn the stopcock in the patient line ON to the patient to obtain pressure and confirm system function.
 8. Zero balance the system regularly per hospital schedule. Check the monitor calibration per hospital protocol.

Elastic Patient Mounting Strap

Use the Elastic Patient Mounting Strap to mount transducer to the arm of the patient. When using an elastic patient mounting strap observe the following:

1. To assemble, insert the housing of the transducer and hook surface of the strap upright. Press the pointed end of the strap through the top of one slot in the side of the transducer. Pull through until the "D" ring is within 3 inches of the transducer. Pass the strap under the transducer and out the slot in the other side. Do not twist the strap.
2. To mount to the patient, position the transducer/strap assembly on the arm. Loop the untwisted strap around the arm and pass the end through the "D" ring. Pull the strap through the ring with enough tension to hold the transducer in place. Do not over-tighten. Loop the strap back over itself and secure with the loop fastener.

Continuous Flush Device

This product will continuously infuse a pressurized solution. The 3 mL/hr device is for adult continuous infusion. It nominally flows at 3 mL/hr with an inlet pressure of 300 mmHg. The 30 mL/hr device (yellow band) is only for use with an infusion pump or other flow controller. It delivers a nominal flow of 30 mL/hr with an inlet pressure of 300 mmHg. When using a continuous flush device, observe the following:

1. Connect the infusion fluid source to the female luer lock connector.
2. To prime the system or supply a fast flow of infusion fluid to the patient, either pull the flush tab or squeeze the wings on the continuous flush device. Release the pull tab or wings to return to the continuous flush mode.

Stopcocks, Manifolds, Monitoring Lines, and Arterial Extensions

When using a stopcock, manifold or arterial extension, observe the following:

1. Inspect components for damage (including kinks).
2. To close a single port, position the stopcock that the word OFF is over the port.
3. To close all ports on a 4-way stopcock, position the handle 45 degrees from either side of any port.
4. Align all stopcock handles exactly on the center of the port for the best possible response.

IV Accessories

The micro-drip chamber nominally gives 60 drops per mL. This produces 3 drops per minute with a 3 mL/hr flow rate. The macro-drip chamber nominal drop size is 20 drops per mL. This produces 1 drop per minute with a 3 mL/hr flow rate.

Needle-less Sample Site

When using the needle-less sampling site observe the following:

1. Use aseptic technique. To access valve port: Swab valve port for 3 seconds and allow to dry prior to use or per institution protocol.
2. Rotate handle to control fluid flow.
3. Flow will be through ports aligned with indicator ribs on handle.
4. Port aligned with handle will be closed.
5. Use all ports open position of the handle in order to enable a continuous flow through the entire stopcock and the luer activated port. Continuous flow occurs from the patient distal port through the valve to the patient proximal port.
6. During IV application use all ports open position of the handle, except for when side port (luer activated valve) is in use.
7. Prior to detaching a syringe from the luer activated valve, handle position should be closed. Make sure syringe is depressurized in order to avoid fluid drop from syringe.
8. Flush the line according to hospital protocol and/or CDC guidelines.

NOTE:

- These instructions are NOT meant to define or suggest any medical or surgical technique. The individual physician is responsible for the proper procedure and techniques to be used with this product.
- When priming the needle sample site make sure handle is at all open position. No male luer should be inserted into the valve during priming.

Hyperbaric Pass Through and Infusion Kits

The Hyperbaric Disposable Pass Through and Intravenous Infusion Extension Kits are intended for use with the Sechrist Monoplace Hyperbaric chamber and appropriate infusion pump. These kits are designed to interface with the chamber bulkhead door facilitating the transfer of IV fluids into a pressurized hyperbaric chamber. When using a hyperbaric pass through and infusion kit observe the following:

1. Inspect the package and seal before opening.
2. Open the package and inspect the product.
3. Grasp the pass through and insert the female end through the desired fitting located on the inside wall of the chamber door. Insert the pass through until fully inserted. Tighten knurled nut on the exterior of the chamber door.
4. If provided, attach stopcock from tray to female port of pass through. Attach other components as necessary and finger tighten.
5. Fluid fill kit. Determine if using normally closed or normally open back check valve. Normally closed back check valve is clear and the normally opened back check valve is white plastic.

Normally Closed Back Check Valve

If using normally closed back check valve, fluid fill the kit using infusion pump or 10 mL syringe.

- With stopcock attached to female port of pass through, attach syringe to side port of stopcock. Attach desired IV source to other port of stopcock. Draw enough IV solution into syringe to fill all branches of kit (6 mL minimum for standard Infusion Kit). Close stopcock to the IV source and inject solution into kit. (Use tray to collect excess fluid.) Turn handle OFF to side port of stopcock. Remove syringe from stopcock. Cap side port of stopcock with plug.
- When not using stopcock, attach desired IV source to female port of pass through and fluid fill using pressure provided by infusion pump. The normally closed back check valve requires a minimum of one psi fluid pressure to activate. System will not gravity fill.

Normally Open Back Check Valve

If using normally open back check valve, attach desired IV source to hyperbaric pass through or port of stopcock. Connect IV extension set to patient's catheter without allowing air into system.

KeepSafe™ Closed Sampling System

KeepSafe™ is a feature added to pressure monitoring kits that is designed to facilitate blood sample collection with reduced risks of external contamination exposure to the patient and blood exposure to the clinician.

1. Swab the surface of the needle-less sample site with the preferred antiseptic for 3 seconds and allow to dry prior to use or per institutional protocol.

2. To clear saline from the needle-less sample site, turn the reservoir stopcock OFF to the transducer. Pull gently on the syringe plunger until the blood reaches the Sampling Site. Draw sufficient blood past the Sampling Site to ensure the sample's purity. This may require adding an extra reservoir syringe to some kits.
3. Turn the handle on the needle-less sample site OFF to the transducer.
4. Draw blood sample with a needle-free syringe through the luer activated side port.
5. Prior to detaching a syringe from the luer activated valve, handle position should be off to the luer activated side port. Make sure syringe is depressurized in order to avoid a drop from the syringe.
6. Use all ports open position of the handle in order to enable a continuous flow through the entire stopcock and the luer activated valve.
7. Remove any residual blood or saline from the syringe reservoir and turn the stopcock OFF to the reservoir.
8. Flush the line according to hospital protocol and/or CDC guidelines.
9. During pressure readings, turn handle to align with the luer activated side port to close this port and prevent potential dampening of pressure readings by the elastomeric valve.

Argo Bagz Waste Management System

1. Connect spike(s) to the intended fluid management system and/or contrast management system.
2. Prime the system with intended fluid.
3. Aspirate intended fluid through the one-way relief valve.
4. Inject waste through the one-way relief valve into waste bag.

Disposal: After use the product may be potential biohazard. Handle and dispose of in accordance with accepted medical practice and all applicable laws and regulations.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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