


A New, CLEANER Indication for Treating Thrombus in the Peripheral Vasculature



The CLEANER™ Rotational Thrombectomy System is now indicated for **mechanical declotting** and controlled and selective **infusion** of physician-specified fluids, including thrombolytics, in the **peripheral vasculature**.

Designed by:

REX MEDICAL®

ARGON
MEDICAL DEVICES



CLEANER XT ORDERING INFORMATION


QUANTITY	CATALOG NUMBER	CATHETER OD	LENGTH	SINUSOIDAL WIRE AMPLITUDE
_____	700009XT	6F	65cm	9mm
_____	700109XT	6F	135cm	9mm



CLEANER 15 ORDERING INFORMATION

QUANTITY	CATALOG NUMBER	CATHETER OD	LENGTH	SINUSOIDAL WIRE AMPLITUDE
_____	700015	7F	65cm	15mm
_____	700115	7F	135cm	15mm

PHYSICIAN'S SIGNATURE



1445 Flat Creek Road
Athens, Texas 75751 USA
800.927.4669 Customer Service
903.677.9396 Fax
www.argonmedical.com

© 2014 Argon Medical Devices, Inc., All Rights Reserved.

No reproduction of any portion of this work shall be made without prior permission of the author.

* CLEANERXT and CLEANER15 are trademarks of Rex Medical

Designed By:



96-6000-36/1114A